

Multi Purpose Unit.

1. Goel Offer multipurpose pilot plant for chemical and pharmaceutical industries for process development, scale-up, process simulation and kiloscale cGMP production in batch and semi-batch operation. The pilot plant used for chemical processing includes solid charging, liquid charging, reaction, heating / cooling, rectification, auto / manual reflux arrangement,

2. The multipurpose pilot plant designed in such a way that we can modify the same easily as per process requirement.

Pilot Plant Available with:

1. Jacketed full glass reactor/ Cylindrical full glass reactor with Oil heating cooling bath / Spherical full glass reactor with Oil heating cooling bath
2. Multipurpose glass distillation overhead
3. Stainless steel / MS epoxy coated / MS painted frame supporting
4. Flame proof / Non flame proof / cGMP / non GMP models available
5. Excellent corrosion resistant.
6. Temp. Controller.
7. Gas purging, solid charging / multi liquid addition.
8. Vacuum / exhaust piping
9. Additional feeders / receivers
10. Solid feeding

Unit Cat. Ref.	Reduction Capacity	Bath KW	Addition Vessel	Vapour Line	Condenser HTA m2	Cooler HTA m2	Receiver Size
MPU 20	20 L	4.5	2 L	80 DN	0.35	0.10	2L, 5L
MPU 50	50 L	6.0	5 L	100 DN	0.50	0.20	5L, 10L
MPU 100	100 L	9.0	10 L	150 DN	1.50	0.35	10L, 20L
MPU 200	200 L	12.0	20 L	150 DL	1.50	0.35	10L, 20L
MPU 300	300 L	16.0	20 L	225 DL	2.50	0.50	20L, 20L

