


JACKETTED VESSELS (TRIPLE WALL)

Goel is presenting Flexi (Detachable) “Transparent Double Jacketed Vessel” to maintain the leading position in industry by way of developing indigenously newer products with higher value towards their Mission and to cater the customer via innovation.

Salient Features of Jacketed Vessel

Vacuum Jacket ensures

- -Transparent insulation.
- -Minimum heat loss.
- -Process visibility.
- -Minimize frost formation for cryogenics operation.
- **Vacuum Jacket & main vessels are detachable ensures**
- -Ease of cleaning.
- -No breakages due to thermal expansion.
- -Partial Replacement of any part can be possible reduce the maintenance cost.



Cat.Ref.	L	L1	L2	A	B	D1	D2	DN	DN1	DN2
VZT5/6/8	620	350	175	75	310	215	275	150	25	25
VZT10/9/12	700	420	175	90	370	315	390	225	25	25
VZT20/12/14	675	500	175	125	368	365	423	300	25	40
VZT50/16/18	825	650	175	200	350	465	523	400	25	40